


MSB

CROWN BRIDGE




Crown Instruments

Kroneninstrumente
Pinze per corone
Alicates para coronas
Pinces pour couronnes




**KERLING
DP-38A**
145 mm

Kupferringabnehmezange
Copper Ring Removing Pliers
Pinza per togliere anelli di rame
Pinza para sacar anillos de cobre
Pince à enlever les bagues en cuivre




**FURRER
DP-38**
145 mm




DP-38R
Schraube
Replacement Screw
Vite di ricambio
Tornillos de recambio
Vis de rechange


Teleskopkronenzangen - Pliers for telescope crowns
Pinze per corone telescopiche - Alicates para coronas telescópicas
Pinces pour couronnes telescopes



**BÖHM
DP-28**
150 mm



**BÖHM
DP-29**
150 mm



DP-24
120 mm



DP-24C
Schlüssel
Key
Chiavetta
Llave
Clé



DP-24R
Diamantierte Einsätze, paar
Diamond inserts, pairs
Punte diamantate, coppia
Puntas diamantadas, par
Pointes diamantées, pair


Crown Instruments

Kronen Instrumente
Strumenti per corone
Instrumentos para coronas
Instruments pour couronnes

Crown and Bridge

Krone und Brücke
Couronne et pont
Corona e Ponte
Korona i most


MRT-C7
Complete Bridge Extractor
MRT-7 + C-71 + C-72 + C-73 + C-74




Removing pliers
for temporary
bridges and crowns




Tungsten-Carbide Inserts
Mit Hartmetalleinsätzen
Plaquettes en carbure de tungstène
Plaquetas de carburo de tungsteno
Inserti in carburo di tungsteno




MRT/C7



MRT-7



MRT-6
170 mm



PH-60
170 mm

Mit Silikon
With silicone insert
Con inserto in silicone
Con inserto de silicone
Avec insert en silicon

1/1

Crown Removers

Kronenentferner


Leva corone

Removedor de coronas


Enlève-couronne

Enlève-couronne


MRT-C3
Complete with 3 Attachments
C-34 + C-35 + C-36



MRT-C4
Complete with 3 Attachments
C-41 + C-42 + C-43




MRT-C1
Complete with 3 Attachments
C-11 + C-12 + C-13



C-30

Usando la punta C-30
e invertendo la molla
si usa come spingicorone



MRT-C3



MRT-C4




MRT-C1

Crown Removers

Kronenentferner
Leva corone
Removedor de coronas
Enlève-couronne

Crown and Bridge

Krone und Brücke
Couronne et pont
Corona e Ponte
Korona i most



UNIVERSAL
DL-60
150 mm




DL-C
140 mm



DL-V
130 mm



DL-0
130 mm



DP-136
Crown Pliers
140mm

TEMPORARY CROWN PLACEMENT KIT | DP-136

1 Crown Pliers + 20 Replacement Tips + 10g Adhesive Powder

Crown Pliers: designed to remove upper and lower crowns and bridges. An adjustable stopper allows to adjust the grasping tips to the width of the crown and locking into place. It will prevent any unwanted pressures to the crown avoiding potential damages. The double spring back action along with the hook-lock mechanism allows for a smooth grasping motion.

Polymer Plastic Replacement Tips: provide a perfect and safely grip for temporary ceramic crowns without damaging metal, resin, or ceramic surfaces. A better grip of the crown can be achieved if the tips are coated with the supplied adhesive powder which is dampened by saliva.

Adhesive Powder: denture adhesive used as "glue" for recementing crowns or bridges that have come off.

DP-137 / Replacement Kit for DP-136: 20 Replacement Tips + 10g Adhesive Powder